

Usage instructions LUXURY Egg Decorating Kit

• Crayon:

1. Boil the eggs and let them cool and dry. The eggs should be warm.
2. Draw on the undyed eggshell with the crayon to create a drawing or a pattern.
3. Dye the eggs following the instructions below. The dye will not color the drawing or the pattern on the egg.

• Egg dye tablets:

1. Put one tablet in a cup.
2. Add 5 oz of warm water in each cup and 2-3 tablespoons of vinegar. Stir and wait until the tablet dissolves.
3. Bend the metal holder and dip the egg.
4. Remove the egg from the solution after obtaining the desired color.
5. Take out the egg and let it dry on the drying tray.

• Stickers:

Use already colored, cool and dry eggs. Gently peel the stickers from the sticker sheet and press them tightly to the eggshell.

• Thermal stickers:

1. Choose medium sized eggs, boil them and leave them to cool. Separate the thermal stickers apart following the dotted line.
2. Place the egg into the thermal sticker so that the tips are seen equally from both ends.
3. Place the egg on a tablespoon or hold it with suitable tongs. Put the egg in boiling water. The egg is ready and tightly wrapped in the thermal sticker in 1-2 seconds.
4. If you want the egg ends to be colored dye the egg after putting the sticker on. Dry well and remove the dye left on the sticker. Thermal stickers are suitable for cracked eggs.

• Egg stands:

Use scissors to cut the egg stands from the sides of the box and assemble them by inserting cutting 1 into cutting 2.

• Drying tray:

Punch out the pre-cut circles from the back of the box and use it as a drying tray.





Egg Decorating Kit **luxury**

Warning and Precautions: Choking hazard. May contain small parts.
Not suitable for children under 3 years of age.
This product may stain some surfaces. Use warm water and soap to clean.

Colors



- YELLOW • ORANGE • RED •
- LIGHT GREEN • GREEN •
- LIGHT BLUE • BLUE •
- PURPLE • PINK •

Contents



41 stickers




12 egg stands



6 thermal stickers



Suitable for white and brown eggs.
For coloring up to 90 eggs.



2 metal holders



9 egg dye tablets



1 crayon

Instructions



Put one tablet in a cup



Add 5 oz of warm water in each cup and 2-3 tablespoons of vinegar. Stir and wait until the tablet dissolves



Bend the metal holder and dip the egg



Remove the egg from the solution after obtaining the desired color



Take out the egg and let it dry

Usage instructions inside!

Use warm and boiled eggs.

Egg Dye Tablets

Ingredients:

Lactose, Microcrystalline Cellulose, FD&C Yellow No.5, FD&C Red No.40, FD&C Red No.3, FD&C Blue No.1, FD&C Blue No.2, Cornstarch, Magnesium Stearate, Silicon Dioxide.

Contains: Milk

Conforms to ASTM D-4236.

Food safe dyes and all other components are made and packed in the European Union.